

Mini



TS3M35

Regular



TS3S40



TS3S45



TS3S50

S-L-A Implant

L

- 7mm
- 8.5mm
- 10mm
- 11.5mm
- 13mm
- 15mm

- 311010
- 311020
- 311030
- 311040
- 311050

- 312010
- 312020
- 312030
- 312040
- 312050
- 312060

- 313010
- 313020
- 313030
- 313040
- 313050
- 313060

- 314020
- 314030
- 314040
- 314050
- 314060
- 314070

- 7mm
- 8.5mm
- 10mm
- 11.5mm
- 13mm
- 15mm

Lance Drill



Φ2.0
351030



Φ2.0
351030



Φ2.0
351030



Φ2.0
351030

Twist Drill

L

- 7mm
- 8.5mm
- 10mm
- 11.5mm
- 13mm
- 15mm

- Φ2.0
351080
- 351110
- 351140
- 351170
- 351250
- 351380

- Φ2.0
351080
- 351110
- 351140
- 351170
- 351250
- 351380

- Φ2.0
351080
- 351110
- 351140
- 351170
- 351250
- 351380

- Φ2.0
351080
- 351110
- 351140
- 351170
- 351250
- 351380

- 7mm
- 8.5mm
- 10mm
- 11.5mm
- 13mm
- 15mm

Direct Drill

L

- 7mm
- 8.5mm
- 10mm
- 11.5mm
- 13mm
- 15mm



	Φ3.0(Soft Bone)	F3.5(Normal Bone)
7mm	351180/351300	351700
8.5mm		351710
10mm		351720
11.5mm		351730
13mm	351300	351740
15mm		351750



	F3.5(Soft Bone)	F4.0(Normal Bone)
7mm	351700	351760
8.5mm		351770
10mm		351780
11.5mm		351790
13mm	351740	351800
15mm		351810



	F3.5	F4.0(Soft Bone)
7mm	351700	351760
8.5mm		351770
10mm		351780
11.5mm		351790
13mm	351740	351800
15mm		351810



	F3.5	F4.5(Soft Bone)
7mm	351700	351820
8.5mm		351830
10mm		351840
11.5mm		351850
13mm	351740	351860
15mm		351870



	F4.5(Normal Bone)
7mm	351820
8.5mm	351830
10mm	351840
11.5mm	351850
13mm	351860
15mm	351870



	F5.0(Normal Bone)
7mm	351880
8.5mm	351890
10mm	351900
11.5mm	351910
13mm	351920
15mm	351930

Taper Cortical Drill

Use after forming a final straight drill hole in normal and hard bone.
Lower line is used in normal bone case.
Upper line is used in hard bone case.



F3.5(Hard Bone)
351610



F4.0(Hard Bone)
351620




F4.5(Hard Bone)
351630



F5.0(Hard Bone)
351640

1.Cover Screw

2.Healing Abutment



Cover Screw	Φ4.0	Φ4.5	Φ5.5
321010(H0.4)	323010(H3.0)	323110(H3.0)	323310(H3.0)
321020(H2.0)	323020(H4.0)	323120(H4.0)	323340(H4.0)
	323030(H5.0)	323130(H5.0)	323320(H5.0)
	323040(H7.0)	323140(H7.0)	323330(H7.0)



Cover Screw	Φ4.0	Φ5.0
322010(H0.4)	324010(H3.0)	324110(H3.0)
322020(H2.0)	324020(H4.0)	324120(H4.0)
	324030(H5.0)	324130(H5.0)
	324040(H7.0)	324140(H7.0)



Cover Screw	Φ5.0	Φ5.5	Φ6.0
324110(H3.0)	324510(H3.0)	324210(H3.0)	324310(H3.0)
324120(H4.0)	324520(H4.0)	324220(H4.0)	324320(H4.0)
324130(H5.0)	324530(H5.0)	324230(H5.0)	324330(H5.0)
324140(H7.0)	324540(H7.0)	324240(H7.0)	324340(H7.0)



Cover Screw	Φ6.0	Φ7.0	Φ7.5
324210(H3.0)	324610(H3.0)	324710(H3.0)	324810(H3.0)
324220(H4.0)	324620(H4.0)	324720(H4.0)	324820(H4.0)
324230(H5.0)	324630(H5.0)	324730(H5.0)	324830(H5.0)
324240(H7.0)	324640(H7.0)	324740(H7.0)	324840(H7.0)

Mini

Regular

Impression Component

Prosthetic Restoration

Abutment

Temporary Abutment /Esthetic-low Healing Cap

Plastic Coping

Auxiliary Parts Screw

Fixture Pick-up Impression Coping

L	D	Φ4.0
11		361110
15		361120
Open tray		

Fixture Transfer Impression Coping

L	D	Φ4.0
11		361160
14		361170
Closed tray		

Esthetic-low Pick-up Impression Coping

Φ4.8

362470

Open tray

Esthetic-low Transfer Impression Coping

Φ4.8

362480

Closed tray

Fixture Pick-up Impression Coping

L	D	Φ4.0	Φ5.0	Φ6.0
11		362210	362220	362230
15		362240	362250	362260
Open tray				

Fixture Transfer Impression Coping

L	D	Φ4.0	Φ5.0	Φ6.0
11		362300	362320	362340
14		362310	362330	362350
Closed tray				

Esthetic-low Pick-up Impression Coping

Φ4.8

362470

Open tray

Esthetic-low Transfer Impression Coping

Φ4.8

362480

Closed tray

Fixture Lab Analog



Esthetic-low Lab Analog



Fixture Lab Analog



Esthetic-low Lab Analog



Cement Retained Restoration

Cement Retained Restoration

Screw Retained Restoration

Cement Retained Restoration

Cement Retained Restoration

Screw Retained Restoration


Transfer Abutment



H	GH	D	Φ4.5
5.5	1.0	Hex	331010
			331020
			331030
			331040
			331050

* Tightening torque: 30Ncm

Angled Abutment



H	GH	Type	D	Φ4.5	
8.0	2.0	A	333010	333020	
					B
		4.0	A	333040	333050
			B	333040	333050

* Tightening torque: 30Ncm


Multi Abutment



GH	D	Φ4.8
1.0	3.0	337010
		337020
		337030
		337040
		337050

* Tightening torque: 30Ncm


Multi-Angled abutment



A	GH	D	Φ4.8
17°	2.5	3.0	337130
			337140
			337150
30°	4.0	3.5	337160
			337170
			337180

* Tightening torque: 25Ncm


Transfer Abutment



H	GH	D	Φ4.5	Φ5.0	Φ6.0	Φ7.0	
5.5	1.0	Hex	Hex	331060	331160	331310	331410
			Hex	331070	331170	331320	331420
			Hex	331080	331180	331330	331430
			Hex	331090	331190	331340	331440
			Hex	331100	331200	331350	331450
7.0	1.0	-	Hex	-	331210	331360	-
			Hex	-	331220	331370	-
			Hex	-	331230	331380	-
			Hex	-	331240	331390	-
			Hex	-	331250	331400	-

* Tightening torque: 35Ncm

Angled Abutment



H	GH	Type	D	Φ5.0	Φ6.0		
8.0	2.0	A	30°	Hex	333070	333130	
				Hex	333080	333140	
		4.0	B	30°	Hex	333100	333160
					Hex	333110	333170

* Tightening torque: 35Ncm


Multi Abutment



GH	D	Φ4.8
1.0	3.0	337070
		337080
		337090
		337100
		337110

* Tightening torque: 30Ncm


Multi-Angled abutment



A	GH	D	Φ4.8
17°	2.5	3.0	337190
			337200
			337210
30°	4.0	3.5	337230
			337240
			337250

* Tightening torque: 30Ncm

Temporary Abutment



H	GH	Type	D	Φ4.0	
10	1.0	Hex	371010	371020	
					Non-hex
		3.0	Hex	371030	371040
				Non-hex	371040

* Tightening torque: 20-30Ncm

Esthetic-low Temporary Cylinder




Esthetic-low Healing Cap



* Tightening torque: 20Ncm

Temporary Abutment



H	GH	Type	D	Φ4.5	
10	1.0	Hex	371050	371060	
					Non-hex
		3.0	Hex	371070	371080
				Non-hex	371080

* Tightening torque: 20-35Ncm

Esthetic-low Temporary Cylinder



Esthetic-low Healing Cap



* Tightening torque: 20Ncm

Esthetic-low Plastic Cylinder



Esthetic-low Cylinder Screw



377020

* Tightening torque: 20Ncm

Esthetic-low Plastic Cylinder

Type	Crown	Bridge
Φ4.5	371310	371360
Φ5.0	371320	371370
Φ6.0	371330	371380
Φ7.0	371340	371390

347010

* Tightening torque: 20Ncm

Esthetic-low Plastic Cylinder



Esthetic-low Cylinder Screw



377020

* Tightening torque: 20Ncm

Lab Screw



Waxing Screw



Lab Screw




Waxing Screw




Finishing Reamer Set



Hex Hand Drill




Hex Torque Drill



Multi Abutment Outer Driver




Multi-angled Abutment Carrier




Finishing Reamer Set




Hex Hand Drill



Hex Torque Drill



Multi Abutment Outer Driver



Multi-angled Abutment Carrier



Esthetic-low Polishing Protector



Esthetic-low Polishing Protector

